

Protein Interactions in the Mycofactocin Biosynthetic Pathway

Protein-protein interactions are essential for many biochemical transformations and cellular processes. However, traditional methods of detecting these interactions are often times qualitative and provide little informationⁱ. Alternatively, more robust quantitative methods require large quantities of valuable protein sample if the dissociation constant is highⁱⁱ or rely on innate protein fluorescence or fluorescence labeling for detectionⁱⁱⁱ. Surface plasmon resonance (SPR) is capable of analyzing protein-protein interactions kinetically while consuming little protein sample and does not rely on absorbance/fluorescence properties of proteins. This application note describes the SPR analysis of protein-protein interactions found in the mycofactocin biosynthetic pathway.

The radical-S-methionine (RS) protein MftC belongs to a subfamily of proteins known to modify peptides and through gene knockout studies, it was found to be critical for *M. tuberculosis* growth with cholesterol as the sole carbon source^{iv}. Recently, it was shown that peptide chaperones, or small proteins of ~100 amino acids in length, play an essential, yet not fully understood role in this process^v. In the first step of mycofactocin biosynthesis, the peptide chaperone MftB was found to be critical for MftC to catalyze the modification of the peptide MftA. To measure the dissociation constants between MftB and MftC, a five channel SPR (BI-4500) and a nickel-nitrilotriacetic acid (Ni-NTA) Au sensor chip was used. As shown in Figure 1, the process involves activating the NTA chip with Ni²⁺, followed by binding the His₆-tagged protein to the Ni-NTA chip, exposing the analyte to the ligand, and stripping Ni²⁺ and the protein complex with ethylenediaminetetraacetic acid (EDTA) to regenerate the bare NTA chip.




Figure 1 A schematic representation of an assay using His₆-tagged protein and a Ni-NTA Au Chip. The assay involves activating the chip with NiCl₂, binding the his-tagged protein, exposing the analyte to the ligand, and stripping the Ni-protein complex with EDTA.

In the assay, His₆-MftC were bounded to three channels and exposed to various concentrations (0.5, 1.0, 2.5, 5.0, and 10 μM) of MftB at a flow rate of 60 μL/min (Figure 2A). Background was subtracted from all data using a channel exposed to MftB without MftC present. Analysis of the steady state parameters (Figure 2B), where the maximum change in response units (RU) is plotted against the concentration of MftB, yielded a K_D = 2.2 ± 0.3 μM^{vi}. An independent kinetic analysis was used to verify the steady-state measured dissociation constant. Three concentrations (1, 2.5, and 5 μM) were chosen and kinetically analyzed each trace individually (Figure 2C) to determine their k_a and k_d rates and thereby generating the K_D

$(K_D = k_d/k_a)^{vi}$. The kinetic analysis yielded a $K_D = 1.3 \pm 0.7$, in good agreement with the steady-state analysis. In summary, by using the SPR technique, it is shown that protein-protein interactions are present in the mycofactocin biosynthetic pathway.

