

Binding studies of an Integrin and its Mutants using BI-4500

Integrins are heterodimeric receptors which help in cell adhesion and are ubiquitous to all metazoans.¹ They also play a vital role in numerous cellular processes including cell growth, proliferation, differentiation, migration, and phagocytosis of opsonized target.^{2,3} Integrin, Mac-1 is expressed primarily on myeloid cells such as neutrophils and macrophages. It is involved in the phagocytosis of opsonized particles, cell-cell fusion, and acts as an important regulator of the immune activities of these cells.⁴ To date, more than 100 ligands have been reported to bind to Mac-1.²

The binding site for most ligands is in the I-domain of Mac-1's α M subunit (~ 200 residue α MI-domain). Despite the importance of α MI-domain in ligand binding, only the structure of α MI-domain with one of its natural ligands, the complement fragment iC3b, has been determined thus far.^{5,6} Similar to other integrins, Mac-1 can adopt both inactive and active conformations. The active conformation has an increasing affinity for its ligands; most studies of α MI-domain's interactions with ligands use active α MI-domain.⁷ Ligand specificities of integrins are vital determinants of their activities, and understanding the ligand binding mechanisms behind these specificities is critical to develop drugs targeting integrins.




Figure 1: Schematic of SPR measurement of Mac-1 receptor α MI-domain.

In this study, three mutants of α MI-domain were purified and characterized using SPR⁸ as seen in Figure 1. Among the three mutants (Δ K315, I316G, and Q163C/Q309C) only one of the mutants (Q163C/Q309C) can spontaneously form the required disulfide bond to maintain the active conformation of the α MI-domain and has a better yield with thermal stability. Studies claim that HEK293T cells expressing Mac-1 with the Q163C/Q309C mutations were shown to have a higher affinity for ligands than wild type Mac-1 expressing cells, but there is no exclusive biochemical study showing the Q163C/Q309C α MI-domain has enhanced affinity for ligands compared to wild type I-domain.⁹ Surface Plasmon Resonance (SPR) was used to understand the affinity of the isolated mutants with C3d and to investigate the enhanced affinity of Q163C/Q309C mutants.

SPR analysis was carried out on a BI-4500A SPR instrument (Biosensing Instrument Inc.). Figure 2 shows the SPR sensorgrams for the α MI-domain variants flowed over a sensor with immobilized C3d. These sensorgrams fitted to a one-to-one binding model showed that, Q163C/Q309C, Δ K315, and I316G mutants all have dissociation constant of binding (Kd) like other reported studies. The results indicated that the Q163C/Q309C mutant has a similar affinity for C3d as other active mutants. However, the koff rate of Δ K315 appears to be slower than the other two active mutants, implying the kinetics of binding differ slightly among the mutants. Control experiments carried out using inactive α MI-domain showed inactive α MI-domain does not have a measurable affinity for C3d.




Figure 2: SPR characterization of C3d's interactions with active α MI-domains.

Before measurement, 50 μ M C3d was flowed at a rate of 20 μ L/min over EDC/NHS-activated CM-dextran sensor until a response of \sim 1000 RU was observed. An increase of 1 RU corresponds to protein deposition of \sim 1 pg / mm². A) Sensorgrams of a C3d-functionalized sensor treated with 0.19, 0.39, 0.78, 1.56, 3.13, 6.25, and 12.5 μ M of Δ K315, I316G, or Q163C/Q309C mutants. The sensor was regenerated with 1.5 M NaCl after each sample. Injection of α MI-domain started at time zero and ended after 180 seconds. B) Estimation of the dissociation constant (K_d) of the interaction for each active α MI-domain by fitting the equilibrium values of the response curves to a one-to-one binding model. Response curves were background corrected by subtracting the response of the reference channel with no immobilized C3d.

The integrin Mac-1 is an important integrin involved in many aspects of leukocyte biology. Understanding the mechanisms of its ligand specificity is essential to developing targeted treatments against Mac-1. However, the lack of a suitable active α MI-domain has so far prevented investigating the interactions of active α MI-domain with its ligands using BI-4500A SPR instrument. This report systematically examined known active mutants of α MI-domain and showed that the Q163C/Q309C mutant adopts the active conformation and has an enhanced affinity for ligands compared to wild type I-domain. The availability of such a mutant will enable more SPR studies of α MI-domain's interaction with ligands and reveal more insights into the mechanisms of Mac-1 activity.

The integrin Mac-1 is an important integrin involved in many aspects of leukocyte biology. Understanding the mechanisms of its ligand specificity is essential to developing targeted treatments against Mac-1. However, the lack of a suitable active α MI-domain has so far prevented investigating the interactions of active α MI-domain with its ligands using BI-4500A SPR instrument. This report systematically examined known active mutants of α MI-domain and showed that the Q163C/Q309C mutant adopts the active conformation and has an enhanced affinity for ligands compared to wild type I-domain. The availability of such a mutant will enable more SPR studies of α MI-domain's interaction with ligands and reveal more insights into the mechanisms of Mac-1 activity.

References:

- 1) Luo BH et al; *annurev. immunol.*25.022106.141618 PMID: 17201681
- 2) Lamers C et al; PMID: 33995387
- 3) Vorup-Jensen T et al; PMID: 30534123
- 4) Rosetti F et al; PMID: 26683153
- 5) Bajic G et al; *pnas.*1311261110 PMID: 24065820
- 6) Fernandez FJ et al; *Nat Commun* 13: 1955.
- 7) Lee JO et al; PMID: 8747460
- 8) Nguyen et al; *Plos one* 18.1 (2023): e0280778.
- 9) Shimaoka M et al; *PNAS* 99: 16737–16741.