

## Binding studies of an Integrin and its Mutants using BI-4500

Integrins are heterodimeric receptors which help in cell adhesion and are ubiquitous to all metazoans.<sup>1</sup> They also play a vital role in numerous cellular processes including cell growth, proliferation, differentiation, migration and phagocytosis of opsonized target.<sup>2,3</sup> Integrin, Mac-1 is expressed primarily on myeloid cells such as neutrophils and macrophages. It is involved in the phagocytosis of opsonized particles, cell-cell fusion, and acts as an important regulator of the immune activities of these cells.<sup>4</sup> To date, more than 100 ligands have been reported to bind to Mac-1.<sup>2</sup>

The binding site for the vast majority of ligands is located in the I-domain of Mac-1's  $\alpha$ M subunit (~ 200 residue  $\alpha$ MI-domain). Despite the importance of  $\alpha$ MI-domain in ligand binding, only the structure of  $\alpha$ MI-domain with one of its natural ligands, the complement fragment iC3b, has been determined thus far.<sup>5,6</sup> Similar to other integrins, Mac-1 can adopt both inactive and active conformations. The active conformation has an increasing affinity for its ligands; most studies of  $\alpha$ MI-domain's interactions with ligands use active  $\alpha$ MI-domain.<sup>7</sup> Ligand specificities of integrins are vital determinants of their activities, and understanding the ligand binding mechanisms behind these specificities is critical to develop drugs targeting integrins.





Figure 1: Schematic of SPR measurement of Mac-1 receptor  $\alpha$ MI-domain.

In this study, three mutants of  $\alpha$ MI-domain were purified and characterized using SPR<sup>8</sup> as seen in Figure 1. Among the three mutants ( $\Delta$ K315, I316G, and Q163C/Q309C) only one of the mutants (Q163C/Q309C) can spontaneously form the required disulfide bond to maintain the active conformation of the  $\alpha$ MI-domain and has a better yield with thermal stability. Studies claim that HEK293T cells expressing Mac-1 with the Q163C/Q309C mutations were shown to have a higher affinity for ligands than wild type Mac-1 expressing cells, but there is no exclusive biochemical study showing the Q163C/Q309C  $\alpha$ MI-domain has enhanced affinity for ligands compared to wild type I-domain.<sup>9</sup> Surface Plasmon Resonance (SPR) was used to understand the affinity of the isolated mutants with C3d and to investigate the enhanced affinity of Q163C/Q309C mutants.

SPR analysis was carried out on a BI-4500A SPR instrument (Biosensing Instrument Inc.). Figure 2 shows the SPR sensorgrams for the  $\alpha$ MI-domain variants flowed over a sensor with immobilized C3d. These sensorgrams fitted to a one-to-one binding model showed that, Q163C/Q309C,  $\Delta$ K315, and I316G mutants all have dissociation constant of binding ( $K_d$ ) similar to other reported studies. The results indicated that the Q163C/Q309C mutant has a similar affinity for C3d as other active mutants. However, the  $k_{off}$  rate of  $\Delta$ K315 appears to be slower than the other two active mutants, implying the kinetics of binding differ slightly among the

mutants. Control experiments carried out using inactive  $\alpha$ MI-domain showed inactive  $\alpha$ MI-domain does not have a measurable affinity for C3d.



**Figure 2: SPR characterization of C3d's interactions with active  $\alpha$ MI-domains.** Before measurement, 50  $\mu$ M C3d was flowed at a rate of 20  $\mu$ L/min over EDC/NHS-activated CM-dextran sensor until a response of  $\sim$  1000 RU was observed. An increase of 1 RU corresponds to protein deposition of  $\sim$  1 pg / mm<sup>2</sup>. A) Sensorgrams of a C3d-functionalized sensor treated with 0.19, 0.39 0.78, 1.56, 3.13, 6.25, and 12.5  $\mu$ M of  $\Delta$ K315, I316G, or Q163C/Q309C mutants. The sensor was regenerated with 1.5 M NaCl after each sample. Injection of  $\alpha$ MI-domain started at time zero and ended after 180 seconds. B) Estimation of the dissociation constant ( $K_d$ ) of the interaction for each active  $\alpha$ MI-domain by fitting the equilibrium values of the response curves to a one-to-one binding model. Response curves were background corrected by subtracting the response of the reference channel with no immobilized C3d.

The integrin Mac-1 is an important integrin involved in many aspects of leukocyte biology. Understanding the mechanisms of its ligand specificity is essential to developing targeted treatments against Mac-1. However, the lack of a suitable active  $\alpha$ MI-domain has so far prevented investigating the interactions of active  $\alpha$ MI-domain with its ligands using BI-4500A SPR instrument. This report systematically examined known active mutants of  $\alpha$ MI-domain and showed that the Q163C/Q309C mutant adopts the active conformation and has an enhanced affinity for ligands compared to wild type I-domain. The availability of such a mutant will enable more SPR studies of  $\alpha$ MI-domain's interaction with ligands and reveal more insights into the mechanisms of Mac-1 activity.

## References:

- 1) Luo BH, Carman CV, Springer TA (2007) Structural basis of integrin regulation and signaling. *Annu Rev Immunol* 25: 619–647. <https://doi.org/10.1146/annurev.immunol.25.022106.141618> PMID: 17201681
- 2) Lamers C, Pluss CJ, Ricklin D (2021) The Promiscuous Profile of Complement Receptor 3 in Ligand Binding, Immune Modulation, and Pathophysiology. *Front Immunol* 12: 662164. <https://doi.org/10.3389/fimmu.2021.662164> PMID: 33995387
- 3) Vorup-Jensen T, Jensen RK (2018) Structural Immunology of Complement Receptors 3 and 4. *Front Immunol* 9: 2716. <https://doi.org/10.3389/fimmu.2018.02716> PMID: 30534123
- 4) Rosetti F, Mayadas TN (2016) The many faces of Mac-1 in autoimmune disease. *Immunol Rev* 269: 175–193. <https://doi.org/10.1111/imr.12373> PMID: 26683153
- 5) Bajic G, Yatime L, Sim RB, Vorup-Jensen T, Andersen GR (2013) Structural insight on the recognition of surface-bound opsonins by the integrin I domain of complement receptor *Proc Natl Acad Sci USA* 110: 16426–16431. <https://doi.org/10.1073/pnas.1311261110> PMID: 24065820
- 6) Fernandez FJ, Santos-Lopez J, Martinez-Barricarte R, Querol-Garcia J, Martin-Merinerio H, et al. (2022) The crystal structure of iC3b-CR3 alpha reveals a modular recognition of the main opsonin iC3b by the CR3 integrin receptor. *Nat Commun* 13: 1955.
- 7) Lee JO, Bankston LA, Arnaout MA, Liddington RC (1995) Two conformations of the integrin A-domain (I-domain): a pathway for activation? *Structure* 3: 1333–1340. [https://doi.org/10.1016/s0969-2126\(01\)00271-4](https://doi.org/10.1016/s0969-2126(01)00271-4) PMID: 8747460
- 8) Nguyen, Hoa, Tianwei Jing, and Xu Wang. "The Q163C/Q309C mutant of  $\alpha$ MI-domain is an active variant suitable for NMR characterization." *Plos one* 18.1 (2023): e0280778.
- 9) Shimaoka M, Lu C, Salas A, Xiao T, Takagi J, et al. (2002) Stabilizing the integrin  $\alpha$ M inserted domain in alternative conformations with a range of engineered disulfide bonds. *Proceedings of the National Academy of Sciences* 99: 16737–16741.