

Surface Plasmon Resonance Analysis of Chemically Synthesized Cyclic Malaria Antigen PfAMA1 Domain I Binding to PfRON2

Plasmodium falciparum is responsible for most severe malaria cases, causing over half a million deaths annually. Central to erythrocyte invasion is the interaction between apical membrane antigen 1 (*PfAMA1*) and rhoptry neck protein 2 (*PfRON2*).¹ The *PfRON2* (~4 kDa) ectodomain docks into a conserved hydrophobic cleft on *PfAMA1* domain I (DI, ~30 kDa), and because no alternative invasion pathway exists, disrupting this interaction is an attractive antimalarial strategy. Chemical synthesis of *PfAMA1* offers unique advantages over recombinant expression, including incorporation of non-natural amino acids and critical access to the D-enantiomeric mirror-image protein for biological display screening of D-peptide inhibitors.²

In this study, truncated 180-residue *PfAMA1*-DI construct (His123–Cys302) was chemically synthesized via a convergent six-segment native chemical ligation strategy. An AffiTag—comprising polyhistidine and biotin connected through a flexible 16-residue polyArg-containing spacer was appended to the cyclic construct to enable gentle Ni-NTA surface capture for SPR.³ SPR experiments were performed on the 5-channel BI-4500 instrument.




Figure 1: SPR sensorgrams from three independent replicates (A, B, C) of AffiTag-cyclicPfAMA1-DI binding to PfRON22021–2059 at 600 nM, 800 nM, and 1.0 μ M. Solid orange lines represent 1:1 Langmuir fits to each concentration series. The average KD across the three replicates was 500.24 ± 14 nM.

Initial attempts to immobilize cyclicPfAMA1-DI directly on a dextran sensor chip via amine coupling, and separately to capture biotinylated PfRON2 on a streptavidin (SA) chip, both failed to yield reliable binding data — the former showing no interpretable signal and the latter exhibiting substantial non-specific binding of cyclicPfAMA1-DI to the streptavidin surface itself, likely due to disruption of the protein's functional fold or surface-driven artefacts. These limitations motivated the switch to the gentler, orientation-controlled Ni-NTA/His-tag capture strategy. Folded AffiTag-cyclicPfAMA1-DI was captured on a nickel-NTA sensor chip via its polyhistidine tag; *PfRON2* was passed as analyte at different concentrations (500 nM, 800 nM, and 1.0 μ M). Binding was concentration-dependent and saturable (**Figure 1**). Global fitting to a 1:1 Langmuir model (Figure 2) yielded an average association rate constant $k_a = (1.42 \pm 0.14) \times 10^5$ $M^{-1}s^{-1}$, average dissociation rate constant $k_d = (7.08 \pm 0.81) \times 10^{-2}$ s^{-1} , and average equilibrium dissociation constant $K_D = 500 \pm 14$ nM (**Figure 1, Table 1**).

Measurement run	k_a ($M^{-1}s^{-1}$)	k_d (s^{-1})	K_D (nM)
A	1.43×10^5	6.95×10^{-2}	486.10
B	1.26×10^5	6.34×10^{-2}	501.16
C	1.55×10^5	7.95×10^{-2}	513.48
Average	$(1.42 \pm 0.14) \times 10^5$	$(7.08 \pm 0.81) \times 10^{-2}$	500.24 ± 14

Table 1: Kinetic parameters for each of the three independent SPR replicates (A, B, C) of PfRON22021–2059 binding to AffiTag-cyclicPfAMA1-DI, with the resulting average $\pm$ SD ($n = 3$).

The entire binding experiment was independently repeated three times (**Figure 1**, panels A–C), and each replicate was fitted separately to a 1:1 Langmuir model to extract k_a , k_d , and K_D (Table 1); the three independent estimates were then averaged to give the final reported values. Across the three replicates, k_a ranged from 1.26×10^5 to $1.55 \times 10^5 \text{ M}^{-1}\text{s}^{-1}$ and k_d ranged from 6.34×10^{-2} to $7.95 \times 10^{-2} \text{ s}^{-1}$, giving individual K_D estimates of 486, 501, and 513 nM which supports the reproducibility of the measurement.




Figure 2: Structural comparison of the cyclic PfAMA1 Domain I (DI) construct and the PfAMA1 Domain I + Domain II (DI+DII) model. (a) Predicted structure of the cyclic PfAMA1-DI bound to PfRON2. (b) Overlay of the cyclic PfAMA1-DI model with the PfAMA1 DI+DII structure showing the absence of Domain II (gray) in the cyclic DI construct. Domain I is shown in green, Domain II in gray, and PfRON2 in red/blue.

PfAMA1-DI construct kinetic parameters were compared to the recombinant full-length PfAMA1-(DI+DII) (**Figure 2**). The comparison of kinetic parameters (**Table 1**) reveals that the k_a of cyclic PfAMA1-DI is essentially identical to PfAMA1-(DI+DII) — 1.42×10^5 vs. $1.49 \times 10^5 \text{ M}^{-1}\text{s}^{-1}$ respectively, demonstrating that encounter-complex formation is driven entirely by domain I and is independent of domain II. By contrast, the dissociation rate is ~22-fold faster in the truncated construct (7.08×10^{-2} vs. $3.21 \times 10^{-3} \text{ s}^{-1}$), accounting for the ~23-fold weaker ($K_D=22 \text{ nM}$) of 500 nM. This kinetic penalty is attributable to the absence of domain II, particularly the DII loop (Ala346–Lys395), which has been shown to prolong complex half-life 18-fold by kinetically locking the bound PfRON2 conformation through transient contacts with its N-terminal α -helix. Also, the K_D of 500 nM also aligns closely with the affinity of a truncated PfRON2 variant lacking the N-terminal helix (~520 nM), providing independent structural corroboration.⁴

Despite the weaker absolute affinity relative to the full-length protein, cyclicPfAMA1-DI retains a functional PfRON2-binding groove and represents a viable scaffold for synthesis of the D-enantiomeric mirror-image protein, enabling mirror-image biological display campaigns to identify D-peptide or L-protein inhibitors of the PfAMA1–PfRON2 interaction. The Ni-NTA SPR format described here provides a reproducible, orientation-controlled kinetic assay for such drug discovery programs.

Materials is a summary of Mannuthodikayil et al, Chemistry–A European Journal 31, no. 28 (2025): e202500894

References:

- 1) Tonkin, M. L et al, Science **2011**, 333 (6041), 463–467.
- 2) Kent, S. B. H et al, Chem. Soc. Rev. **2009**, 38 (2), 338–351.
- 3) Mannuthodikayil et al, Chemistry–A European Journal 31, no. 28 (2025): e202500894
- 4) Biswas et al, Biochemistry and Biophysics Reports 26 (2021): 100950.