

Increase your SPR research efficiency with the BI Continuous Flow Syringe Pump. The BI syringe pump PP-100 provides continuous flow during experiments for improved productivity and unattended operation. With its automatic in-line degasser, your SPR system now becomes a fully automated system.

Key Features:

- ✧ Automated in-line degassing
- ✧ Automatic buffer exchange up to six sources
- ✧ Fully integrated with BI software with operational wizards
- ✧ Supplies buffer for entire SPR and Autosampler systems
- ✧ For use with BI-4500A SPR and SPRm 200A



Specifications

Fluid handling	Buffer selection	Automatic exchange between 6 sources
	In-line degassing	Included
	Flow rate	300 uL/min
System	Outer dimension	300 (w) x 190 (h) x 360 (d) mm
	Weight	6.8 kg
	Power	100 - 240V, 50-60Hz, 5A max

Email: Info2BI@BiosensingUSA.com
 Website: www.BiosensingUSA.com
 Tel: 1-480-491-2777 Fax: 1-866-897-8741